



# CYMEDICA E-VIVE™ SIZING GUIDE

<b>e-vive NMES System w/o Brace:</b>				
Conductive Garment, NMES Controller, USB Charger, Electrodes & Electrode Gel				
Garment size	Quad Circumference 8" (20 cm) up from top of patella		Left	Right
	inches	cm		
Extra Small	14.0" - 18.0"	36 - 46	CY-1000-055	CY-1000-054
Small/Medium	18.0" - 22.0"	46 - 56	CY-1000-116	CY-1000-113
Large/X-Large	22.0" - 26.0"	56 - 66	CY-1000-117	CY-1000-114
2X-Large/3X-Large	26.0" - 30.0"	66 - 76	CY-1000-118	CY-1000-115
<b>e-vive NMES System with Brace:</b>				
Post Operative Brace, Conductive Garment, NMES Controller, USB Charger, Electrodes & Electrode Gel				
Garment size	Quad Circumference 8" (20 cm) up from top of patella		Left	Right
	inches	cm		
Extra Small	14.0" - 18.0"	36 - 46	CY-1000-048	CY-1000-047
Small/Medium	18.0" - 22.0"	46 - 56	CY-1000-110	CY-1000-107
Large/X-Large	22.0" - 26.0"	56 - 66	CY-1000-111	CY-1000-108
2X-Large/3X-Large	26.0" - 30.0"	66 - 76	CY-1000-112	CY-1000-109
<b>e-vive Bilateral NMES System:</b>				
2 Conductive Garments, NMES Controller, USB Charger, Electrodes & Electrode Gel				
Garment size	Quad Circumference 8" (20 cm) up from top of patella		System includes both Left and Right Garments	
	inches	cm		
Extra Small	14.0" - 18.0"	36 - 46	CY-1000-050	
Small/Medium	18.0" - 22.0"	46 - 56	CY-1000-051	
Large/X-Large	22.0" - 26.0"	56 - 66	CY-1000-052	
2X-Large/3X-Large	26.0" - 30.0"	66 - 76	CY-1000-053	

If patient is of athletic build and measures near limit or between sizes, it is recommended to select the larger size. Conversely, if patient is in between sizes with poor muscle tone a smaller size may be most appropriate.

[www.cymedicaortho.com](http://www.cymedicaortho.com)  
 19120 N. Pima Rd, Ste. 135  
 Scottsdale, AZ 85255  
 844.296.2014