



CYMEDICA E-VIVE™ SIZING GUIDE

e-vive NMES System w/o Brace:				
Conductive Garment, NMES Controller, USB Charger, Electrodes & Electrode Gel				
Garment size	Quad Circumference 8" (20 cm) up from top of patella		Left	Right
	inches	cm		
Extra Small	15.0" - 18.0"	38 - 46	CY-1000-222	CY-1000-221
Small/Medium	18.0" - 22.0"	46 - 56	CY-1000-216	CY-1000-213
Large/X-Large	22.0" - 26.0"	56 - 66	CY-1000-217	CY-1000-214
2X-Large/3X-Large	26.0" - 30.0"	66 - 76	CY-1000-218	CY-1000-215
e-vive NMES System with Brace:				
Post Operative Brace, Conductive Garment, NMES Controller, USB Charger, Electrodes & Electrode Gel				
Garment size	Quad Circumference 8" (20 cm) up from top of patella		Left	Right
	inches	cm		
Extra Small	15.0" - 18.0"	38 - 46	CY-1000-220	CY-1000-219
Small/Medium	18.0" - 22.0"	46 - 56	CY-1000-210	CY-1000-207
Large/X-Large	22.0" - 26.0"	56 - 66	CY-1000-211	CY-1000-208
2X-Large/3X-Large	26.0" - 30.0"	66 - 76	CY-1000-212	CY-1000-209
e-vive Bilateral NMES System:				
2 Conductive Garments, NMES Controller, USB Charger, Electrodes & Electrode Gel				
Garment size	Quad Circumference 8" (20 cm) up from top of patella		System includes both Left and Right Garments	
	inches	cm		
Extra Small	15.0" - 18.0"	38 - 46	CY-1000-250	
Small/Medium	18.0" - 22.0"	46 - 56	CY-1000-251	
Large/X-Large	22.0" - 26.0"	56 - 66	CY-1000-252	
2X-Large/3X-Large	26.0" - 30.0"	66 - 76	CY-1000-253	

If patient is of athletic build and measures near limit or between sizes, it is recommended to select the larger size. Conversely, if patient is in between sizes with poor muscle tone a smaller size may be most appropriate.

www.cymedicaortho.com
 19120 N. Pima Rd, Ste. 135
 Scottsdale, AZ 85255
 844.296.2014