



# CyMedica e-vive™ Sizing Guide



<b>e-vive NMES System w/o Brace:</b>				
Conductive Garment, NMES Controller, USB Charger, Electrodes & Electrode Gel				
Garment size	Quad Circumference 8" (20 cm) up from top of patella		Left	Right
	inches	cm		
Small/Medium	18.0"- 22.0"	46 - 56	CY-1000-316	CY-1000-313
Large/X-Large	22.0"- 26.0"	56 - 66	CY-1000-317	CY-1000-314
2X-Large/3X-Large	26.0"- 30.0"	66 - 76	CY-1000-318	CY-1000-315
<b>e-vive NMES System with Brace:</b>				
Post Operative Brace, Conductive Garment, NMES Controller, USB Charger, Electrodes & Electrode Gel				
Garment size	Quad Circumference 8" (20 cm) up from top of patella		Left	Right
	inches	cm		
Small/Medium	18.0"- 22.0"	46 - 56	CY-1000-310	CY-1000-307
Large/X-Large	22.0"- 26.0"	56 - 66	CY-1000-311	CY-1000-308
2X-Large/3X-Large	26.0"- 30.0"	66 - 76	CY-1000-312	CY-1000-309
<b>e-vive Bilateral NMES System:</b>				
2 Conductive Garments, NMES Controller, USB Charger, Electrodes & Electrode Gel				
Garment size	Quad Circumference 8" (20 cm) up from top of patella		System includes both Left and Right Garments	
	inches	cm		
Small/Medium	18.0"- 22.0"	46 - 56	CY-1000-351	
Large/X-Large	22.0"- 26.0"	56 - 66	CY-1000-352	
2X-Large/3X-Large	26.0"- 30.0"	66 - 76	CY-1000-353	

If patient is of athletic build and measures near limit or between sizes, it is recommended to select the larger size. Conversely, if patient is in between sizes with poor muscle tone a smaller size may be most appropriate.

[www.cymedicaortho.com](http://www.cymedicaortho.com)  
 19120 N. Pima Rd, Ste. 135  
 Scottsdale, AZ 85255  
 844.296.2014